



## Visual-Visual


You have selected the Visual-Visual task.




Press **'Next'** to receive further instructions.

## *Step 1: Straighten Head*


- If the halo is red you must reorient your head:



- Look straight ahead so that the virtual laser pointer falls on the central target, then slowly tilt your head on your shoulders. This will make the halo turn progressively from red to green:



- When the halo pops to bright green, hold this position.



## *Step 2: Acquire the Target Orientation*

- Once your head is aligned with your body (i.e. when the halo is bright green), a line of balls will appear at the end of the tunnel, representing the target orientation.



- Memorise the target orientation.

### *Step 3: Tilt the Head*

- If the halo turns red, you must tilt your head to a new orientation. The arrow tells you which way to rotate the head.




- Roll your head from side to side until the halo turns bright green and hold this position.



### ***Step 4: Align to Remembered Target***

- With the left hand use the ***VR Headset Remote*** to align the virtual tool to the memorised target.

**Note:** Press and hold the right and left arrows on the ***VR Headset Remote*** to rotate the virtual tool clockwise and counter-clockwise, respectively.



- Press the ***Select Button*** on the ***VR Headset Remote*** to validate your response once you think that the virtual tool is aligned with the memorised target.

### *Step 5: Check Response*

- In some of the trials a line of target balls will be projected from your hand towards the targets, allowing you to verify whether or not you correctly aligned your hand.

